
Rexton Bluetooth Wireless Device Manual

*Rexton
Bluetooth
Wireless
Device
Manual*

*Downloaded from
db.whlarchitecture.com
by Hardy & Fernanda*

INTRODUCTION

For those who rely on hearing assistance technology, the **Rexton Bluetooth Wireless Device Manual** is an essential resource that unlocks the full potential of your hearing aids and wireless accessories. Rexton has established itself as a trusted brand in audiology, offering devices that bridge the gap between clear hearing and modern

connectivity. Whether you have just received your first Rexton hearing aid or you are looking to pair your device with a smartphone or TV streamer, understanding the setup and operation of your Bluetooth wireless device is crucial. This comprehensive guide serves as an extended companion to the official manual, providing step-by-step instructions, troubleshooting tips, and best practices to help you get the most out of your Rexton

Bluetooth wireless device. We will walk you through everything from initial pairing to advanced features like app integration and maintenance, ensuring your experience is seamless and satisfying.

UNDERSTANDING YOUR REXTON BLUETOOTH WIRELESS DEVICE

Before diving into the manual, it is helpful to recognize the components that make up a Rexton Bluetooth wireless system. Typically, this includes a pair of rechargeable or battery-powered hearing aids equipped with Bluetooth Low Energy (BLE) technology, along with optional accessories like the Rexton Remote Microphone+, the TV Streamer+, or the My

Rexton companion app on your smartphone. The Rexton Bluetooth Wireless Device Manual usually covers the basic operation of the hearing aids themselves and how they connect wirelessly to external devices. The core of the system is the ability to stream audio directly from your phone, tablet, or television, allowing you to hear calls, music, and TV sound with remarkable clarity.

It is important to note that different models may have slightly different pairing procedures. However, the foundational steps remain consistent across the Rexton product lineup. The manual typically highlights a simple button press or the use of the My Rexton app to initiate pairing.

Always check your specific model number—often found on the charging case or inside the battery compartment—to ensure you are following the correct instructions. The **Rexton Bluetooth Wireless Device Manual** is designed to be user-friendly, but having additional context can smooth out any rough edges during setup.

GETTING STARTED: UNBOXING AND INITIAL SETUP

What's in the Box?

When you first open your Rexton Bluetooth device package, you will typically find the hearing aids, a charging case, a USB

charging cable, a user manual, and possibly several ear tips or domes. Some packages also include a cleaning tool and a wireless remote or streamer. It is a good idea to inspect all items and verify that nothing is missing. The **Rexton Bluetooth Wireless Device Manual** will list all included components and their functions.

Charging First

Before attempting to pair your Bluetooth devices, ensure your hearing aids are fully charged. If you have a rechargeable model, place them in the charging case and connect the case to a power source using the provided USB cable.

Most Rexton models require about 2-3 hours for a full charge. A blinking LED on the case indicates that charging is in progress. Once the lights are steady green, your devices are ready. For battery-powered aids, insert fresh batteries and close the battery compartment. The manual emphasizes that initial charging or battery installation should be done before any Bluetooth pairing.

Turning On and Basic Operation

To turn on your Rexton hearing aids, simply remove them from the charging case if they

are rechargeable, or close the battery compartment if they are battery-operated. You will hear a start-up melody. The **Rexton Bluetooth Wireless Device Manual** often describes how to adjust volume and change programs using the tactile buttons on the hearing aids or the app. For first-time users, it is advisable to start with a low volume and gradually increase it in a quiet environment.

PAIRING YOUR REXTON BLUETOOTH DEVICE WITH A SMARTPHONE

Step-by-

Step Pairing Process

Pairing your Rexton hearing aids with a smartphone is straightforward, thanks to Bluetooth. The **Rexton Bluetooth Wireless Device Manual** outlines the following general steps:

1. **Download the My Rexton**

App: Go to the App Store (iOS) or Google Play Store (Android) and install the official My Rexton app. This app is your control center for volume, programs, and Bluetooth

settings.

2. **Enable Bluetooth on Your Phone:**

Make sure your phone's Bluetooth is turned on and that it is discoverable.

3. **Initiate Pairing Mode on the Hearing Aids:**

Depending on your model, this may involve pressing and holding the button on the hearing aid or the charging case for a few seconds. Some models also enter pairing mode automatically when you first turn them on. Look for a blinking blue or white LED as

- confirmation.
4. **Select the Device in Your Phone's Bluetooth Menu:** Under available devices, you should see an entry like "Rexton [model name]". Tap on it to pair. You may be asked to confirm a PIN code - usually 0000 or 1234 - but many Rexton devices use Just Pair technology that skips this step.

5. **Confirm in the My Rexton App:** Open the app and follow the on-screen instructions to complete the setup. The app will usually recognize the

hearing aids and guide you through volume adjustments and program settings.

If you encounter any issues during pairing, the manual suggests resetting the Bluetooth connection by turning off both the hearing aids and the phone, then restarting the process from step one. In many cases, removing the hearing aids from the phone's Bluetooth memory and re-adding them solves connection problems.

Pairing with Multiple

USING THE MY Devices

Rexton Bluetooth devices typically support multipoint pairing, meaning you can connect to both your phone and a TV streamer simultaneously. The **Rexton Bluetooth Wireless Device Manual** explains that once a primary device is paired, you can add a secondary device by repeating the pairing process while the primary is connected. However, audio will only stream from one source at a time. For example, if you are listening to TV and a call comes in, the hearing aids will automatically switch to your phone call and then return to TV when the call ends.

REXTON APP

Key Features of the App

The My Rexton app is a powerful tool that extends the functionality of your Bluetooth wireless device. According to the **Rexton Bluetooth Wireless Device Manual**, the app allows you to:

- **Adjust Volume:**
Use the on-screen slider to increase or decrease the sound level in each ear independently.
- **Change**

Programs:

Switch between listening programs such as “Quiet,” “Noise,” “Music,” and “TV” to adapt to different environments.

- **Check Battery**

Status: See the remaining charge for each hearing aid and the charging case.

- **Locate My**

Hearing Aids: If you misplace your devices, the app can help you find them by playing a tone (if within range).

- **Stream Audio:**

The app acts as a gateway for streaming music, podcasts, and phone calls directly to the hearing aids.

- **Customize**

Settings:

You can create personalized profiles for different locations and save them for quick access.

The app also provides access to the digital manual and firmware updates. Regularly checking for updates is recommended to ensure your devices run smoothly.

Pairing the App with Your Device

After downloading the app, it will ask for permission to access Bluetooth and location services. Grant these

permissions to allow the app to discover your hearing aids. The app's pairing wizard will guide you through the same steps as the manual's pairing instructions. Once paired, the app will automatically connect every time it is opened, as long as the hearing aids are within range (approximately 10 meters).

CONNECTING TO A TV STREAMER OR REMOTE MICROPHONE

TV Streamer + Setup

The Rexton TV Streamer+ is a small device that plugs into your TV's audio output

(optical or analog). It transmits TV sound directly to your hearing aids without disturbing others. The **Rexton Bluetooth Wireless Device Manual** details the following setup:

1. Connect the TV Streamer+ to your TV using the provided cables (optical or 3.5 mm audio cable).
2. Plug the Streamer into a power outlet using the USB power adapter.
3. Put your hearing aids in pairing mode (refer to the manual).
4. Press the pairing button on the TV Streamer+ — typically a small button with a Bluetooth icon. The LED will blink to indicate pairing mode.

5. Once the hearing aids detect the Streamer, the connection will be confirmed with a tone. The LED on the Streamer will become steady.

You can now enjoy TV audio privately. Volume and program changes can be made via the app or the hearing aid buttons. The Streamer also has a charging cradle for your hearing aids, making it an all-in-one solution.

Using the Remote Micropho

ne+

The Rexton Remote Microphone+ is beneficial in noisy environments or for one-on-one conversations. It acts as a wireless microphone that you can clip onto the speaker's lapel. The **Rexton Bluetooth Wireless Device Manual** explains that pairing is similar to the TV Streamer: place the remote microphone into pairing mode and then put your hearing aids into pairing mode. When connected, the remote microphone captures the speaker's voice clearly and streams it to your hearing aids, reducing background noise. It can also be used as a table microphone for small group

discussions.

MAINTENANCE AND CARE FOR YOUR REXTON BLUETOOTH DEVICE

Cleaning and Drying

Hearing aids are sensitive to moisture and earwax. The manual recommends wiping the hearing aids daily with a soft, dry cloth. Use the provided cleaning tool to gently remove debris from the microphone ports and receivers. Never use water or alcohol-based cleaners, as they can damage the electronics. Store your Rexton devices in the charging case when

not in use, as the case provides protection and maintains battery health.

Battery Care

For rechargeable models, avoid exposing the hearing aids to extreme temperatures. The **Rexton Bluetooth Wireless Device Manual** advises charging them at room temperature. Over time, battery capacity may diminish, but most batteries last several years. If you notice significantly reduced runtime, consult your audiologist about battery replacement. For disposable battery models, remove the batteries if you plan to store the hearing aids for an extended period

to prevent leakage.

Firmware Updates

Rexton occasionally releases firmware updates to improve performance and connectivity. The My Rexton app will notify you when an update is available. Ensure your hearing aids are fully charged and stay within range of your phone during the update process. The **Rexton Bluetooth Wireless Device Manual** emphasizes not to interrupt the update, as it could damage the device.

TROUBLESHOOTING COMMON BLUETOOTH ISSUES

No Sound or Intermittent Sound

If your hearing aids are paired but you hear no sound, first check the volume level on both the app and the hearing aids. Make sure the correct program is selected. If you are streaming from a phone, verify that the media volume is turned up. Restarting both the phone and the hearing aids often resolves intermittent sound. If the issue persists, unpair and re-pair the devices.

Pairing Failures

A common frustration is when the hearing aids do not appear in the phone's Bluetooth list. The **Rexton Bluetooth Wireless Device Manual** suggests the following:

- Ensure the hearing aids are turned on and in pairing mode (blinking LED).
- Turn off Bluetooth on other nearby devices to avoid interference.
- Delete any old "Rexton" entries from your phone's Bluetooth memory.
- Restart your phone and try

again.

If you still cannot pair, your hearing aids may need to be reset. Instructions for a factory reset vary by model, but typically involve holding down the button for 15-20 seconds until the LED flashes rapidly. After resetting, you must go through the full pairing process again.

App Connection Problems

If the My Rexton app fails to connect, ensure that location services are enabled (required by Android for Bluetooth scanning). Close the app

completely and reopen it. If the app still does not detect your hearing aids, try unpairing the Bluetooth connection from your phone's settings and then using the app's pairing wizard to re-establish the link. The manual also recommends clearing the app's cache or reinstalling it as a last resort.

**ADVANCED TIPS
AND BEST
PRACTICES**

Using Rexton Bluetooth with a Smartwatch

Some users have successfully connected their Rexton